

UNIVERSITÄT LEIPZIG

Beerinformatics

February 12rd 2019, Bled

Jörg Fallmann, Stephan Bernhart





ONCE UPON A TIME



UNIVERSITÄT LEIPZIG

UNIVERSITÄT

LEIPZIG

BEERINFORMATICS



UNIVERSITÄT

LEIPZIG

BRING SCIENCE BACK FROM THE STREETS

C StreetSo	cienceCommur	nity / StreetScienceCon	Watch → 3 ★ Star 2 ^V / ₈ Fork							
<> Code	() Issues 0	ן Pull requests 0 🛛 🕅 F	Projects 0 III Insights							
Bring science to the "street" https://streetscience.community										
T 113 commits		<i>پا</i> 2 branches	♥ 0 releases	🚨 6 contributors	a <u>t</u> a MIT					
Branch: mas	ter 👻 New pull re	quest	Create new file Upload files F	Find file Clone or download *						

Beerinformatics |

THE PLAN



UNIVERSITÄT LEIPZIG

THE SPONSOR

GERMAN NETWORK FOR BIOINFORMATICS INFRASTRUCTURE

THE SAVIOUR



INTRODUCTION

SETTING

WIMP analysis Mapper Comparison Replicates and yeast composition

Conclusion and Outlook



MOTIVATION

- BEER
- Use state of the art techniques in practical course to motivate and entertain our students
- ++ BEER

LINIVERSITÄT

I FIP7IG

- Compare nanopore results to published Illumina results BEER
- Get better results than what was achieved in F.....g
- ?? did we mention BEER

BEERINFORMATICS



UNIVERSITÄT LEIPZIG

RESEARCH TOPICS

- Why is segemehl the best mapper for nanopore reads?
- How reliable is nanopore's WIMP software?
- Why is segemehl the best mapper for nanopore reads?
- How do replicates compare, what kind of yeast do we have?
- How do results generated by 30 students in 10 days compare to results generated by experts us (on a train ride)?

INTRODUCTION

SETTING

WIMP analysis Mapper Comparison Replicates and yeast composition

Conclusion and Outlook



EXPERIMENTAL SETUP

- Cooperate with real Biochemists (The Moerl Lab)
- BEER was chosen for expected diversity and local interests
 - Gose from Bayrischer Platz (Leipzig "Beer" variant with lactic acid explicably never really caught on outside Leipzig)
 - Weissbier (Schneider Weisse Tap 7) to keep the bavarian experementalist in chief interested
 - ▶ Waldbier because wild cherry was no longer available it was wild pear
- Extract DNA from BEER samples
- Use Oxford Nanopore sequencing on extracted DNA
- Use students to analyse data

Beerinformatics | How to buy beer

SAMPLING

- Gose and Schneider Weisse were accquired from Getränkefeinkost Leipzig
- Waldbier was ordered via the internetz from the producer at https://kiesbye.at



WALDRIES

DNA EXTRACTION

 DNA etraction was performed at Mario Mörl's Lab at the Institute for Biochemistry



Beerinformatics | Who's fault is it if something went wrong?

SAMPLE PREPARATION

 Sample preparation was performed at Mario Mörl's Lab at the Institute for Biochemistry



SAMPLE PREPARATION

- ... or at the Institute for Bioinformatics



NANOPORE SEQUENCING

- Institute for Biochemistry or Institute for Bioinformatics



NANOPORE SEQUENCING RESULTS



UNIVERSITÄT LEIPZIG

INTRODUCTION

SETTING

WIMP analysis Mapper Comparison Replicates and yeast composition

Conclusion and Outlook



Beerinformatics | nanopore's What's in my pot

WIMP

- On-the-fly analysis of nanopore sequencing
- Uses centrifuge metagenoics program suite
- Result is annotated taxonomic tree of reads
- Question How reliable is it?

Lineage
 Taxa - Not Classified
 Taxa - Classified

Figure: WIMP taxonomic tree for Gose run 1, 3%

UNIVERSITÄT LEIPZIG



Figure: WIMP taxonomic tree for Gose run 1, 1%



Figure: WIMP taxonomic tree for Gose run 1, 0.5%

UNIVERSITÄT LEIPZIG



Figure: WIMP taxonomic tree for Gose run 1, 0.1%

UNIVERSITÄT LEIPZIG



Figure: WIMP taxonomic tree for Gose run 1, all reads



Figure: WIMP taxonomic tree for Gose run 2, 1%



Figure: WIMP taxonomic tree for Gose run 2, 0.5%



Figure: WIMP taxonomic tree for Gose run 2, 0.1%



Figure: WIMP taxonomic tree for Gose run 2, all reads

UNIVERSITÄT LEIPZIG



Figure: WIMP taxonomic tree for Waldbier 1, all reads

UNIVERSITÄT LEIPZIG



Figure: WIMP taxonomic tree for Waldbier 2, all reads



Figure: WIMP taxonomic tree for Tap 7 2, all reads

UNIVERSITÄT LEIPZIG



Figure: Representative example of long branch.

- Reads from homo sapiens are suspicious
- Might be contamination
- Contamination at our side would be easily distinguishable

- Reads from homo sapiens are suspicious
- Might be contamination
- Contamination at our side would be easily distinguishable





UNIVERSITÄT LEIPZIG



UNIVERSITÄT LEIPZIG

------ THE "BLACK BOX" TOOL WIMP

MAPPED SPECIES

organism	Gose1	Gose2	Waldbier1	Waldbier2	Weißbier
Talaromyces	4 (3)	4 (4)	96 (704)	40 (135)	86 (19)
Saccharomyces cerevisiae	2382 (1812)	1492 (1200)	67117 (58778)	21089 (17798)	27376 (29224)
Saccharomyces eubayanus	351 (55)	244 (39)	43520 (38241)	14179 (11031)	9360 (165)
Bacillus subtilis	1 (-)	-	1 (-)	6 (-)	31 (10)
Cellulosimicrobium	52 (48)	19 (17)	1065 (921)	259 (213)	196 (165)
Geobacillus	-	-	1 (-)	1 (-)	2 (-)
E. coli	-	1 (-)	57 (3)	13 (1)	42 (20)

UNIVERSITAT

43

Beerinformatics

TODO

ES BLASTET DER FALLMANN AM RAUSCHIGEN BACH BLAST[®] » blastn sulte » RID-6257XNJU016



MAPPER COMPARISON

- Compare results of state of the art mappers
- Work out why segemehl outperforms the others
- (Optimize parameters)

OS and hardware

Segemehl runs nativ on linux-OS and after a long ancient yakutian shamanic dance on Windows UbuntuShell (tested on Windows 10)

On Windows 10 with 8 threads segemehl used 95% of RAM (15,1 Gb from 15,9 GB) and 58% of CPU (Intel Core i7) and consumed approx. 62 GB of the paging file



UNIVERSITÄT LEIPZIG

Warning

Beware of cats, they can press ctrl+Z while running over the keyboard and destroy hours of work at 99% done!



UNIVERSITÄT LEIPZIG

------ SEARCH FOR MAPPERS | for Nanopore reads

	Sample	Mapping Tools	number of Reads	Time for Mapping	Time for Indexing
	Gose1	segemehl	2565	60 s	6 min
MADDING GTATICTICS		minimap2	2565	14 s	← inclusive
MAPPING STATISTICS		GraphMap	2565	66 s	54 s
		BLAT	2565	35 s	- inclusive
	Gose2	segemehl	1667	60 s	6 min
		minimap2	1667	15 s	 inclusive
		GraphMap	1667	40 s	54 s
		BLAT	1667	31 s	- inclusive
	Mix	segemehl	34301	2 h 58 min	6 min
Elapsed time for		minimap2	34301	157 s	← inclusive
indexing and menuing		GraphMap	34301	42 min	54 s
indexing and mapping		BLAT	34301	8 h 54 min	- inclusive
	Waldbier	segemehl	106885	1 h 35 min	6 min
		minimap2	106885	128 s	← inclusive
		GraphMap	106885	39 min	54 s
		BLAT	106885	13 h 20 min	- inclusive
	Weissbier	segemehl	40498	120 s	6 min
		minimap2	40498	19 s	 inclusive
		GraphMap	40498	19 min	54 s
		BLAT	40498	8 min 48 s	 inclusive

UNIVERSITAT Master of Science Bioinformatics

UNIVERSITÄT

LEIPZIG

17

UNIVERSITÄT

LEIPZIG

STUDENTS' RESULTS



UNIVERSITÄT

LEIPZIG

STUDENTS' RESULTS





44

Beerinformatics | How do replicates compare, what kind of yeast do we have?

Replicates

- We have two "technical" replicates (Waldbier, Gose)
- We treat BEER types as biological replicates
- How do they differ?
- Focus on yeast composition



UNIVERSITÄT LEIPZIG



UNIVERSITÄT LEIPZIG

Taxonomic composition



UNIVERSITÄT LEIPZIG

Beerinformatics | It is (not) complicated..

OUR RESULTS



UNIVERSITÄT LEIPZIG

WHAT DO THEY WANT TO TELL US?



UNIVERSITÄT LEIPZIG

INTRODUCTION

SETTING

WIMP analysis Mapper Comparison Replicates and yeast composition

Conclusion and Outlook



CONCLUSION

- Organoleptic testing of BEER can be dangerous to your health
- Students are more cost efficient than post-docs
- If you want to enjoy your beer, you probably should not analyze it

WHY IS METAGENOMICS SO HARD



UNIVERSITÄT LEIPZIG

WORSE WAYS TO WASTE BEER



Feuerwehr Aschbach Markt

UNIVERSITÄT LEIPZIG

ACKNOWLEDGEMENTS

We want to thank:

LINIVERSITÄT

I FIP7IG

- Vereinigung von Förderern und Freunden der Universität Leipzig e. V.
- The Mario Moerl Lab
- FE, BB, BG, PV, MM and **RB** in F.
- the students of the practical course
- teaching staff of the practical course (Zasha, Jan, Sandra, Alexander, Prof. Sonja Prohaska)
- His Wholiness Prof Peter F.. Stadler







UNIVERSITÄT LEIPZIG

Thank You!

Jörg Fallmann, Stephan Bernhart

Transcriptome Bioinformatics Group Institute for Bioinformatics

berni/fall@bioinf.uni-leipzig.de
http://www.bioinf.uni-leipzig.de/